



GARDskin Dose-Response

Skin sensitisation potency measurement and risk assessment without a requirement for animals

IDEA Workshop 2025


Andy Forreryd, SenzaGen AB

Introduction – Testing for skin sensitization

OECD Test Guidelines are mapped to the AOP

AOP - Adverse Outcome Pathway

NAM - New Approach Methods (KE 1-3)



GARDskin was the first assay included in an OECD Test Guideline to combine genomics and machine learning for addressing a regulatory endpoint

The GARD technology platform – How it works

Transcriptomic read-out of the biological response

Biological system: Dendritic-like cell line (KE3)

Readout: Gene expression (genes and toxicity pathways)

Dendritic-like
cell line

ex: h-CLAT (CD86/CD54)

Cellular
responses



GARDskin

Gene expression of biomarker signatures
GARDskin: 196 genes.

Sensitizer


Non-sensitizer

Full transparency: Identities of genes being measured available in peer-reviewed scientific literature.

See for example: Johansson et al. (2011) A genomic biomarker signature can predict skin sensitizers using a cell-based in vitro alternative to animal tests. BMC Genomics.

The GARD technology platform – How it works


Genes cover mechanistically relevant toxicity pathways




How to GARD[®] your products in 6 Steps

1


GARD Input Finder



Grown SenzaCells




Add different concentrations of the test substance to the cells




Determine the concentration of the test substance where the cells react and 90% survive

2

GARD Main Stimulation



Take test substance at determined concentration and add to fresh batch of cells



Quality control of the cells

3


RNA extraction



Extract RNA from the cells

4

Gene expression profiling



Importantly: All genes contribute to a final classification, but with different weights

Prediction algorithm:

$$DV = b + \sum_{i=1}^n w_i x_i$$

n: number of variables (n for GARDskin:196)

b: constant (SVM intercept)

w_i : weight for variable i

x_i : Normalized gene expression data for variable i

Prediction model:

Mean DV ≥ 0 : Skin sensitiser (UN GHS category 1)


Mean DV < 0 : Non-sensitiser.

The OECD approval of GARDskin

Extending the applicability domain of NAM-based assays in OECD TG 442

OECD Test Guideline No. 442 E - *In Vitro* Skin Sensitisation

KE 3 in the AOP for skin sensitization: DCs activation



Performance statistics:

GARDskin accuracy:	94%
WLR	82.1- 88.9%
BLR	92%

Breaking new ground: first harmonised method integrating genomic data for a regulatory endpoint.

Validation study published in peer-reviewed scientific journal:

GARDskin: Published in Johansson et al. (2019), Validation of the GARD™skin assay for assessment of chemical skin sensitizers - ring trial results of predictive performance and reproducibility. *Toxicological Sciences*.

Hydrophobic compounds

Lubrizol: Forreryd, A et al. (2022). Exploration of the GARD™ skin applicability domain: Indirectly acting haptens, hydrophobic substances and UVCBs.

Fragrances & Fragrance formulations

RIFM/IFF: GARDskin dose-response assay and its application in conducting quantitative risk assessment (QRA) for fragrance materials using a NGRA framework.

Metals and metal salts

Johnson & Matthey: Forreryd, A. et al. (2022). The GARD™skin assay: Investigation of the applicability domain for metals.

Complex mixtures

Corteva: Corvaro, M., et al. (2024). GARD™skin and GARD™potency: a proof-of-concept study to investigate the applicability domain for agrochemical formulations.

Polymeric materials

DSM: *In vitro* assessment of skin sensitizing potential of process-related impurities in polymeric materials during product development. Poster.

UVCBs & Natural Complex Substances

Exxon: Assessing the Utility of the GARDskin Assay to Detect Dermal Sensitization Potential in UVCBs and Formulated Lubricant Products.

Surfactants

Work in progress – to be updated.

GARD

Genomic Allergen Rapid Detection

GARDskin Dose-Response: Quantitative assessment of skin sensitizing potency.




GARDskin Dose-Response

Quantitative assessment of skin sensitizing potency

- Perform the GARDskin assay in a titrated range of concentrations ($n \geq 6$).
- Apply standard GARDskin protocol to generate a decision value (DV) for each concentration.
- Estimate cDV_0 : lowest concentration required to induce a positive classification ($DV \geq 0$).

- The GARDskin cDV_0 value correlates significantly with both human NESIL values and LLNA EC3.
- A simple linear regression model allows for continuous potency predictions of a NESIL-value in dose per unit area ($\mu\text{g}/\text{cm}^2$) and can be used directly as a PoD for risk assessment.




GARDskin Dose-Response

Quantitative assessment of skin sensitizing potency

The GARDskin cDV_0 value correlates significantly with both human NESIL values and LLNA EC3.

- Intrinsic limitations exist in both human NESIL and LLNA reference data.
- Both references inform on the same phenomenon i.e., skin sensitizing potency, but neither is perfect and associated with measurement errors.
- It was considered redundant to fit models separately to LLNA EC3 and human NESIL.

→ A composite potency score integrate information from both LLNA and human data into a single value equivalent to a predicted NESIL (in the unit $\mu\text{g}/\text{cm}^2$).



Example	LLNA EC3 ($\mu\text{g}/\text{cm}^2$)	Human NESIL ($\mu\text{g}/\text{cm}^2$)	Composite ($\mu\text{g}/\text{cm}^2$)
DNCB	13.5	8.8	
Cinnamic aldehyde	250	591	
Citral	1450	1420	
...			


GARDskin Dose-Response

Quantitative assessment of skin sensitizing potency


The GARDskin cDV_0 value correlates significantly with both human NESIL values and LLNA EC3.

- Intrinsic limitations exist in both human NESIL and LLNA reference data.
- Both references inform on the same phenomenon i.e., skin sensitizing potency, but neither is perfect and associated with measurement errors.
- It was considered redundant to fit models separately to LLNA EC3 and human NESIL.

→ A composite potency score integrate information from both LLNA and human data into a single value equivalent to a predicted NESIL (in the unit $\mu\text{g}/\text{cm}^2$).



Example	LLNA EC3 ($\mu\text{g}/\text{cm}^2$)	Human NESIL ($\mu\text{g}/\text{cm}^2$)	Composite ($\mu\text{g}/\text{cm}^2$)
DNCB	13.5	8.8	9.8
Cinnamic aldehyde	250	591	378
Citral	1450	1420	1440
...			



GARDskin Dose-Response


Quantitative assessment of skin sensitizing potency

A simple linear regression model allows for continuous potency predictions of a NESIL-value in dose per unit area ($\mu\text{g}/\text{cm}^2$).

- The correlation between cDV_0 and potency is described by a linear regression model.
- The regression model is simple and only contains 1 parameter:

$$\text{Prediction in } \mu\text{g}/\text{cm}^2 = \text{cDV}_0 \text{ in } \mu\text{g}/\text{ml} \times \theta$$

- The model use cDV_0 as input to predict the skin sensitizing potency of an unknown chemical in dose per unit area ($\mu\text{g}/\text{cm}^2$).



GARDskin Dose-Response

How to derive continuous potency predictions


Step 1: Dose-Response testing.

- Generation of a dose response curve by plotting DV vs concentration.
- Identification of a cDV_0 value using linear interpolation.

Step 2: Continuous potency predictions.

- Correlation between cDV_0 and potency is described by a linear regression model.
- The cDV_0 value is used as input into the regression model to derive a potency prediction in the unit $\mu\text{g}/\text{cm}^2$ (LLNA EC3/Human NESIL)

Test Item: Benzyl Cinnamate



GARDskin Dose-Response

How to derive continuous potency predictions


Step 1: Dose-Response testing.

- Generation of a dose response curve by plotting DV vs concentration.
- Identification of a cDV_0 value using linear interpolation.

Step 2: Continuous potency predictions.

- Correlation between cDV_0 and potency is described by a linear regression model.
- The cDV_0 value is used as input into the regression model to derive a potency prediction in the unit $\mu\text{g}/\text{cm}^2$ (LLNA EC3/Human NESIL)

Test Item: Benzyl Cinnamate



GARDskin Dose-Response

Performance for potency predictions


Background

- The performance and reproducibility of GARDskin DR have been evaluated in collaboration with industry (3 papers, ~200 chemicals).

Results – performance

- GARDskin DR predicted human NESIL values correlated well with human reference data (n=40, r=0.69).
- GARDskin DR predicts human NESILs **equally well or better** than the LLNA (LLNA vs Human: r= 0.6).

Note: It is challenging to find reliable reference data. The performance values above refer to a subset of data for which highly curated reference data from both LLNA and human studies were available, as documented in OECD TG 497, Annex 2.



GARDskin Dose-Response

Performance for potency predictions


Results - reproducibility

- GARDskin DR data available for repeated measurements (n=27).
- The predicted human NESIL values were reproducible, with an expected variation between runs of 1.78-fold-changes.

Conclusion

- GARDskin DR provides accurate and reproducible potency predictions with high correlation to human NESIL values and LLNA EC3.

→ The predicted potency value is equivalent to a human NESIL value ($\mu\text{g}/\text{cm}^2$) and can be used directly as a PoD for QRA (see for example Donthamsetty et al. 2024)



Summary & conclusions

- The GARDskin Dose-Response assay is built on the validated framework of GARDskin (OECD TG 442E).
- GARDskin Dose-Response provides accurate and reproducible potency predictions that can be used to establish a safe concentration of a sensitizer in a formulated product.
- Data from the previous IDEA WS meeting demonstrated high performance for the chemicals in the original RCPL list:
 - Continuous potency predictions from the GARDskin Dose-Response assay correlated well with PVs for chemicals in the RCPL list (pearson correlation: 0.74).
 - Overall, very similar potency rankings with GARDskin Dose-Response and RCPL potency list (spearman: 0.69)

Scan QR code for
access to additional
case studies and
publications



Collaboration L'oreal: Development and refinement of GARDskin Dose-Response.

Collaboration RIFM & IFF: Use in NGRA for risk assessments and to determine the maximum concentration of a sensitizer in various consumer products.

Collaboration Takasago: Use during development of novel fragrance materials to establish a safe dose for confirmatory HRIPT.

Collaboration BIC skin care: Use during development of novel temporary ink tattoos to formulate safe and effective products.

Collaboration DoTerra: Use of GARDskin Dose-Response to predict skin sensitization threshold concentrations for Natural Complex Substances.