

# Analysis of hydroperoxides in consumer products

Progress report on the work under the International Dialogue on the Evaluation of Allergens (IDEA)  
 Andreas Natsch, Susanne Kern, Elise Corbi, Christophe Pérès, Marco Nägelin, Hans Leijts, Michel van Strien,  
 Michael J. Calandra, Ying Wang, André Düsterloh, Ulrika Nilsson, Matthias Vey, Cecile Gonzalez, Elena Giménez-Arnau



## The question:

- High frequency of positive patch test reactions to oxidized terpenes are reported.
- No known consumer exposure source for sensitizing doses of hydroperoxides (HP) from ox. terpenes.
- Fragrance products are under suspicion – **is this assumption correct?**
- This question needs to be resolved by analytical methods

## IDEA analytical taskforce: Mission

1. Provide methods to determine hydroperoxides formed from Linalool and Limonene
2. Validate methods for analysis of consumer products
3. Perform market surveillance to test **consumer/patient exposure** from fragranced products

## 1. Toolbox of methods

- GC-MS reduction method: Reduces HP to alcohols; alcohols are reliably measured in different bases
  - Widely available, robust
  - May lead to some overestimation
- Three LC-based methods: Directly detect parent HP, can be used for confirmatory analysis

## Key research steps:

- **Method evaluation:** methods from 6 laboratories
- Two blind-coded ring trials for pre-selection and evaluation of methods
- **Method validation** by 2 **blinded ring trials**
- **Market surveillance** with validated method on **consumer products and products from patients**


## 2. Validation of methods by blinded ring-trials

- Two studies: a) fine fragrances and b) creams and lotions
- Blind-coded samples
- 4 laboratories used GC-MS reduction method in parallel
- 3 laboratories made confirmatory analysis by different LC-MS methods

**Table 1:** Detection of Linalool-OOH in Eau de Toilette (EdT) by **Toolbox of methods**

(data in µg/ml)	EdT No Spike	EdT Low Spike	EdT High Spike
LC-Q-TOF MS	0.0	90.0	279.0
HPLC-CL	0.0	79.5	310.7
LC-orbitrap-MS	0.2	95.7	398.7
Red-GC-MS	6.1	88.6	325.4
<b>spike level added</b>	<b>0.0</b>	<b>92.0</b>	<b>322.0</b>

**Figure 1:** Reduction method validation in blind-coded ring trial



## 3. Market Surveillance study: Key findings

- 104 products analyzed by independent third party CRO
- Only one product with confirmed HP level > 50 µg/ml (i.e. 90 µg/ml limonene-1-OOH)
- **No HP detected in products from (patch test positive) patients**
- **No evidence for HP accumulation in aged samples**
- All results validated by standard addition

## Market Sample: Dose-per area calculation:

- For the single sample with highest level (91 µg/g), exposure per area calculated
- Compared with reporting limit and with toxicological data and patch test dose

### Example of a patient product

Sample and history of donating patient	Analytical methods	Limonene-1-OOH	Limonene-2-OOH	Linalool-7-OOH	Linalool-6-OOH	
O12, Body cream, Positive some fragrances, Positive Limonene ox	GC-MS red. (µg/ml)	<22	<22	<22	<22	⇐ Reduction method
	GC-MS red. (% recovery)	69%	70%	59%	84%	⇐ Spike recovery
	LC-Orbitrap-MS (µg/ml)	NF	nr	NF	NF	⇐ LC-MS method 1
	LC-Q-ToF-MS (µg/ml)	<5	<5	<5	11	⇐ LC-MS method 2
	LC-CL (µg/ml)	NF	NF	NF	NF	⇐ Chemilum. method

### Dose per area calculations for limonene-1-OOH

	Dose of hydroperoxide in test preparation	Dose per area
LLNA Dose inducing sensitisation (EC3)	3300 µg/g (0.33%)	82.5 µg/cm <sup>2</sup>
Patch test limonene-HP *, routine diagnostic level	3300 µg/g (0.33%)	156 µg/cm <sup>2</sup>
Patch test limonene-1-OOH **, diagnostic level	5000 µg/g (0.5%)	228 µg/cm <sup>2</sup>
<b>Defined reporting limit</b>	50 µg/g	0.1 – 0.5*** µg/cm <sup>2</sup>
<b>Analytical data market surveillance:</b> (Max. value of n = 104)	90 µg/g (0.009%)	0.2 µg/cm <sup>2</sup> ****

\* Mixture of isomers, not specifically 1-OOH-isomer  
 \*\* Dose used in study on specific Limonene-1-OOH isomer by Christensson, Contact Dermatitis 2015  
 \*\*\* Different dose depending on product type (Cream 10 mg/cm<sup>2</sup> higher than fine fragrance, 2.2 mg/cm<sup>2</sup>)  
 \*\*\*\* Based on the typical application dose of fine fragrance per area

## Key conclusions:

- Toolbox of methods established
- GC-MS method validated in blinded ring study
- Market surveillance found very low incidence of positive samples
- Low dose in single confirmed positive sample
- Neither detected level nor incidence of positive samples can explain patch test frequency
- **Exposure source for potential sensitization remains elusive**